

<b>Product Name :</b> Astra Automated Nucleic Acid Purification	<b>Product Code :</b> ASTRALAB-SUPPLYA1015
--	---



<b>Description :</b>
Astra Automated Nucleic Acid Purification
<b>Technical Specification :</b>
There are two types of Auto-Pure series, Auto-Pure16A and Auto-Pure 32A. They are desktop devices. Purify DNA and RNA from different starting materials (such as blood, cultured cells or bacteria, tissues, cell-free body fluids and plant samples) by using magnetic separation technology.
The device can use a magnetic bar to transfer particles in each purification stage of binding, mixing, washing

and elution, thus providing a solution with the shortest hands-on time.

DNA and RNA purification has high quality and high molecular weight.

Auto-pure16A can accept up to 16 samples at the same time, while Auto-Pure 32A can accept up to 20 samples at the same time. The maximum sample size is 1ml.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265