

**Product Name :**  
Astra Food Rapid Test Kit

**Product Code :**  
ASTRALAB-SUPPLYA2004



**Description :**

Astra Food Rapid Test Kit

**Technical Specification :**

The rapid test strip product consists of test strips and related reagents. It is easy to operate, with a test time of 5-10 minutes and a favorable price. It is suitable for rapid and qualitative detection of genetically modified solid components in grain seeds, leaves and various other crops.

The rapid detection kit uses double antibody sandwich enzyme-linked immunosorbent assay to detect

genetically modified proteins in samples, with high sensitivity, stable quality, and a shelf life of up to one year.

Quick test kit specifications: antibody-coated 96-well reaction plate, enzyme label and reagents.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265