

<b>Product Name :</b> Astra Food Rapid Test Kit	<b>Product Code :</b> ASTRALAB-SUPPLYA2004
--	---



<b>Description :</b>
Astra Food Rapid Test Kit
<b>Technical Specification :</b>
<p>The rapid test strip product consists of test strips and related reagents. It is easy to operate, with a test time of 5-10 minutes and a favorable price. It is suitable for rapid and qualitative detection of genetically modified solid components in grain seeds, leaves and various other crops.</p> <p>The rapid detection kit uses double antibody sandwich enzyme-linked immunosorbent assay to detect</p>

genetically modified proteins in samples, with high sensitivity, stable quality, and a shelf life of up to one year.

Quick test kit specifications: antibody-coated 96-well reaction plate, enzyme label and reagents.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265