

**Product Name :**  
Astra Automatic Kjeldahl Analyzer

**Product Code :**  
ASTRALAB-SUPPLYA2009



**Description :**

Astra Automatic Kjeldahl Analyzer

**Technical Specification :**

The automatic nitrogen generator distillation device is based on the classic Kjeldahl method.

It is a special instrument for detecting nitrogen content in seeds, dairy products, beverages, feed, soil and other agricultural and sideline products.

It has a special plastic-sprayed steel shell, which is beautiful and practical; the acid-base contact part is made of mirror stainless steel, which is corrosion resistant and durable

Automatic water replenishment function, the liquid receiving volume can be adjusted arbitrarily

Advanced water-saving function design can save a lot of tap water, and advocate new concepts of energy saving, emission reduction and environmental protection

With the steam generator, the water level is controlled by the all-stainless steel electric heating tube, which is corrosion resistant and durable; the water shortage alarm is automatic

The use of a dedicated integrated circuit control system can meet the safe use of different water quality in different regions, ensure normal operation without changing distilled water for a long time, save a lot of distilled water, reduce test costs, and achieve the goal of fast, accurate, and low consumption

The splash-proof tube is made of special high-performance acid and alkali resistant plastic. Good sealing performance, better alkali resistance than glass, and long service life



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265