

Product Name : Astra Milk Test Kit	Product Code : ASTRALAB-SUPPLYA3005
--	---



Description :
Astra Milk Test Kit
Technical Specification :
It can be used for rapid detection of fat, protein, lactose and non-fat solids. It is the first choice for daily milk analysis laboratory users, covering small-scale and daily detection tasks.
Milk fat content, specific gravity, protein content, lactose, relative density, fat, whole milk solids, DDT, mercury, aflatoxin B1, adulteration test, sugar sensory evaluation, microorganisms, adulterated milk, moisture, solubility

and impurities, etc.

Quick and accurate measurement of milk components;

Low content, no harmful chemicals are used;

Detection speed: the effect can be obtained in a short time.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265