

Product Name :
Astra Small Rotary Evaporator

Product Code :
ASTRALAB-SUPPLYA4004



Description :

Astra Small Rotary Evaporator

Technical Specification :

The small rotary evaporator is a small rotary evaporator with a capacity of 2000ml. The 304 stainless steel water bath is used to provide a heat source with a temperature range of $\pm 2^{\circ}\text{C}$.

Small rotary evaporator can be tilted horizontally in the range of $0-45^{\circ}$.

High evaporation rate of vertical condensation system.

Teflon evaporator tube, fluorine rubber gasket, good sealing performance.

Automatic temperature control digital display.

GG-17 high borosilicate glass evaporating flask and receiving flask.

110 V / 60HZ can be customized.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265