

Product Name :
Astra 10L Rotary Evaporator

Product Code :
ASTRALAB-SUPPLYA4006



Description :

Astra 10L Rotary Evaporator

Technical Specification :

The 10L rotary vacuum evaporator has a double-coil condenser, which can achieve a better condensation effect. It has an explosion-proof version.

Teflon (PTFE) and fluoro rubber double rotary seal, patented technology to ensure high vacuum

The temperature control system is accurate and reliable.

Under vacuum distillation conditions, In does not affect the system and solution. Continuous collection

PTFE discharge valve, corrosion-resistant, no pollution.

Use corrosion and high temperature resistant borosilicate glass.

The rubber parts adopt new standard materials, which is convenient for users to replace.

We have two types: conventional type and explosion-proof type, which can be selected according to the experimental environment.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265