

**Product Name :**  
Astra Hydrothermal Synthesis Reactor

**Product Code :**  
ASTRALAB-SUPPLYA5005



**Description :**

Astra Hydrothermal Synthesis Reactor

**Technical Specification :**

Working temperature: ?230°C

Working pressure: ?2MPa.

It has a complete PTFE material lining and excellent acid and alkali resistance.

The pot body and lid are sealed with a durable and reliable sealing structure, and the sealing effect can be stable for a long time without leakage.

Working temperature: ?230°C

Working pressure: ?2MPa (pressure gauge)

Heating and cooling rate: ?5°C/min

Use high-quality 0Cr18Ni9Ti non-magnetic stainless steel.

The pot body and lid are sealed with a durable and reliable sealing structure, and the sealing effect can be stable for a long time without leakage.

Complete PTFE material lining, with excellent acid and alkali resistance.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265