

**Product Name :**  
Astra Laboratory Hot Air Oven

**Product Code :**  
ASTRALAB-SUPPLYA6003



**Description :**

Astra Laboratory Hot Air Oven

**Technical Specification :**

The laboratory hot air oven has a capacity of 215L. All are made of corrosion-resistant 304 stainless steel. The cast iron cab is under high vacuum pressure without any deformation. The cuboid chamber can maximize the effective volume, and the microcomputer temperature controller ensures accurate and reliable temperature control.

Lock air tightness and high vacuum

The cast iron cab can withstand high vacuum pressure without any deformation

The pressure makes the water boil 7 degrees Celsius

Low humidity ensures drying speed and less time

Energy consumption is 20% lower than the previous model

Cuboid cavity can maximize the effective volume, microcomputer temperature controller, accurate and reliable temperature control.

The studio is made of stainless steel plate (or brushed plate) to ensure durability.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265