

**Product Name :**  
Astra Freeze Dry Machine

**Product Code :**  
ASTRALAB-SUPPLYA6004



**Description :**

Astra Freeze Dry Machine

**Technical Specification :**

The tray spacing of the chemical vertical freeze dry machine is 70mm, which can make the freeze dryer have a lower temperature and a better freeze drying effect. The inside of the freeze dryer equipment is made of 304 stainless steel and has good chemical properties.

Embedded circuit design, 128 M FLASH, can store a large amount of freeze-dried data;

ELarge opening trap, no coil inside, with sample pre-freezing function.

Precision refrigeration compressor, low noise.

The cold trap and control panel are made of stainless steel, which is corrosion resistant and easy to clean.

The inflation valve is configured to be filled with dry inert gas;



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265