

Product Name :
Astra Food Freeze Dryer

Product Code :
ASTRALAB-SUPPLYA6006



Description :

Astra Food Freeze Dryer

Technical Specification :

This product is a vacuum freeze dryer that can be connected to a vacuum pump. Vacuum pumping helps the freeze dryer reach a lower operating temperature. We can provide you with vacuum equipment, such as rotary vane vacuum pumps, desktop circulating water vacuum pumps and diaphragm pumps to meet your needs.

Large trap, no internal coil, no sample pre-freezing

The control system automatically saves freeze-drying data, and can be viewed in real-time and historical curves, and the entire freeze-drying process is very clear.

This machine can store a large amount of freeze-dried data, and use U disk to copy the data to the computer, and use the software to browse and print on the computer.

U disk can store freeze-dried data.

Configure the inflation valve to fill dry inert gas.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265