

Product Name : Astra Food Freeze Dryer	Product Code : ASTRALAB-SUPPLYA6006
--	---



Description :
Astra Food Freeze Dryer
Technical Specification :
<p>This product is a vacuum freeze dryer that can be connected to a vacuum pump. Vacuum pumping helps the freeze dryer reach a lower operating temperature. We can provide you with vacuum equipment, such as rotary vane vacuum pumps, desktop circulating water vacuum pumps and diaphragm pumps to meet your needs.</p> <p>Large trap, no internal coil, no sample pre-freezing</p>

The control system automatically saves freeze-drying data, and can be viewed in real-time and historical curves, and the entire freeze-drying process is very clear.

This machine can store a large amount of freeze-dried data, and use U disk to copy the data to the computer, and use the software to browse and print on the computer.

U disk can store freeze-dried data.

Configure the inflation valve to fill dry inert gas.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265