

**Product Name :**  
Astra Cryogenic Vacuum Pump

**Product Code :**  
ASTRALAB-SUPPLYA7002



**Description :**

Astra Cryogenic Vacuum Pump

**Technical Specification :**

The small chiller is suitable as a condensing device for short path distillation kits and small rotary evaporators. Can help the still work more efficiently.

Cooling chiller has professional refrigeration equipment, protection equipment and capacitors.

Dedicated to refrigeration units, expected relays, protective devices, and capacitors;

It can effectively protect the normal use of all kinds of precision equipment and instruments under the conditions of no water, no water pressure, no water temperature and high temperature working conditions.

The circulatory system is made of anti-corrosion material and has functions such as anti-rust, anti-corrosion, and low-temperature liquid pollution sterilization.

It can be matched with various instruments (rotary evaporators, fermentation tanks, chemical reactors, freeze-drying equipment, pharmaceutical reactors, etc.);

With microcomputer software lock function, it can lock the system parameter settings.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265