

**Product Name :**  
Medical Circulating Water Bath

**Product Code :**  
ASTRALAB-SUPPLYA8001



**Description :**

Medical Circulating Water Bath

**Technical Specification :**

This circulating water bath uses an all-stainless steel housing. Designed and produced a circulating constant temperature heating device. It has the characteristics of compact and reasonable structure, safe and convenient use, precise temperature control, controllable output flow, and energy saving. The volume of circulating water and oil bath is 5L, 10L, 20L, 30L, 50L, 80L, 100L.

The whole machine is made of imported stainless steel plate, with high-quality electrostatic spraying plastic shell.

The electric heating tube is immersed in the oil, heating is fast, the thermal efficiency is high, the electric energy is saved, and the electronic constant temperature.

There is insulation wool between the outer shell and the lining of the instrument, which has excellent heat insulation effect.

The instrument adopts digital display temperature control, simple and eye-catching operation.

PID control, with the advantages of high temperature control accuracy and small temperature marks.

The circulation system is made of stainless steel and copper; it has the characteristics of anti-rust and anti-corrosion.

Solid state relay control circuit, non-contact type, halogen lamp, prolong service life.

Vertical structure with casters for easy movement.

The high-performance circulating pump with circulating water and oil bath provides a constant temperature heat source for the reactor.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265