

Product Name : Medical Humidity Incubator	Product Code : ASTRALAB-SUPPLYA10002
---	--



Description :
Medical Humidity Incubator
Technical Specification :
Internal volume: 250L~2000L;
Temperature range: 10?100?, temperature control fluctuation: $\pm 0.5^{\circ}$, temperature deviation: $\pm 1.0^{\circ}$;
Controller: Microcomputer controller, sensitive in sensing and small system error;

Compressor: Original imported fully enclosed compressor, high efficiency, low noise; to ensure long-term continuous operation of the equipment;

Data printing: It has a data acquisition system that meets GMP requirements, standard stylus printer and U disk storage;

Box material: The inner tank is made of mirror stainless steel 304, the four corners are semi-circular, which is easy to clean, and the shelf space in the box is adjustable;

System protection: Configure an independent over-temperature protection system;

Safety device: Compressor overheating and overpressure overload protection, fan overheating, temperature upper and lower limit deviation alarm;

Working environment: +5?35?

Remote alarm: Standard remote alarm port, optional remote mobile phone SMS alarm (with power failure alarm) (optional);



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265