

<b>Product Name :</b> Medical Cooling Incubator	<b>Product Code :</b> ASTRALAB-SUPPLYA10003
--	--



<b>Description :</b>
Medical Cooling Incubator
<b>Technical Specification :</b>
Cooling incubator,Adopting the latest structural design, using imported high quality and manufacturing technology to ensure long-term continuous operation of the equipment, stable and reliable performance, and suitable for stability testing and storage of temperature-sensitive drugs.
Internal volume: 150L~800L

Temperature control range: 5?, temperature control fluctuation:  $\pm 1^\circ$ , temperature deviation:  $\pm 2.0^\circ$ ;

Temperature controller: Original imported digital display controller, high precision, low drift;

Refrigeration compressor: Original imported fully enclosed compressor, long life, low noise;

Defrosting device: automatic frosting device;

Data printing: A data collection system that complies with GMP requirements, standard stylus printer, and print data can be stored for a long time;

Data storage: U disk storage, which can store temperature history data for more than 1 year;

On-site alarm: On-site temperature deviation sound and light alarm;

Remote alarm: (optional) remote mobile phone SMS alarm (with power failure alarm), multiple devices can share a mobile phone card;

Box material: The inner tank is made of mirror stainless steel 304, and the outer shell is made of high-quality steel plate sprayed with plastic;

Safety devices: compressor overheating protection, compressor overpressure protection, compressor overload protection and fan overheating protection;

Installation power supply: AC 220V $\pm$ 10% 50HZ;

Working environment: +5?35?;



## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265