

<b>Product Name :</b> Medical Dry Bath Incubator	<b>Product Code :</b> ASTRALAB-SUPPLYA10006
---	--



<b>Description :</b>
Medical Dry Bath Incubator
<b>Technical Specification :</b>
The Dry bath incubator uses advanced laser and numerical control processing equipment to produce products; it is used in mining enterprises, laboratories, scientific research units, etc., for drying, baking, wax melting and sterilization.
Internal volume: 70L~1000L;

Temperature control range: RT + 10~200?, constant temperature fluctuation:  $\pm 1.0?$ , temperature accuracy:  $\pm 0.1?$ ;

Temperature controller: digital microcomputer P.I.D temperature controller, with timing function, accurate and reliable temperature control;

Alarm function: independent temperature limit alarm system, which automatically interrupts when the temperature exceeds the limit, and the experiment runs safely without accidents (optional);

Box material: The inner tank is made of mirror stainless steel, and the outer shell is made of high-quality steel plate sprayed plastic;

Timing range: 1-9999min;

Working environment: +5~40?;

Installation power supply: AC 220V $\pm 10?$ 50HZ;



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265