

Product Name :
Medical Cryogenic Freezer

Product Code :
ASTRALAB-SUPPLYA11002



Description :

Medical Cryogenic Freezer

Technical Specification :

High-precision microcomputer temperature control system, which can easily and accurately set the temperature
Platinum resistance temperature sensor can achieve more precise temperature control.

The high-definition and easy-to-operate control panel can set the temperature range from -40°C to -86°C, and

the digital temperature display accuracy is 1°C

Lockable door, keypad lock and password access function can prevent any unauthorized access and operation

Vertical deep freezer, made of galvanized steel inside and painted steel plate outside, and equipped with 4 casters for easy handling

Adjustable 2-3 stainless steel shelves, with inner door and double silicon gasket seal

Chinese patent non-fluorine mixed refrigerant composed of R134a, R404a, R23, R410A, N2, etc.

Use a temperature printer to record 7 days of data. 72-hour backup battery for display of printer and controller

Power supply: 220V /50Hz 1 phase, can be changed as 220V 60HZ or 110V 50/60HZ.

Optional: Chart recorder, CO2 back up system, storage racks/boxes, Remote alarm system



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265