

Product Name : Medical Small Deep Freezer	Product Code : ASTRALAB-SUPPLYA11003
---	--



Description :
Medical Small Deep Freezer
Technical Specification :
The product has complete sound and visual alarms, including high and low temperature alarms, power failure alarms, sensor failure alarms and half-open alarms
This model of refrigerator is equipped with 2~3 high-quality steel wire shelves, which can be adjusted to any height to meet different needs

Standard built-in USB data logger, remote alarm contacts and RS485 interface for monitoring system

There is a cooling fan inside, which can work when the door is closed and stop when the door is open

CFC-free polyurethane foam insulation is environmentally friendly

The electrically heated glass door filled with inserted gas performs well in terms of heat insulation

The medical refrigerator is equipped with 2 sensors. When the main sensor fails, the auxiliary sensor will be activated immediately

The door is equipped with a lock to prevent unauthorized opening and operation



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265