

Product Name : Medical Cooling Chiller	Product Code : ASTRALAB-SUPPLYA11004
--	--



Description :
Medical Cooling Chiller
Technical Specification :
Large capacity Cooling Chiller, the lowest temperature is -120°C. It is a laboratory refrigerator that can handle various occasions.
This cryogenic coolant circulation pump is equipped with expected relays, protective devices and capacitors for refrigeration devices. With super-cooling ability and fast cooling speed, greatly improving work efficiency

Dedicated to refrigeration devices, expected relays, protection devices and capacitors;

Large cooling capacity and fast cooling speed greatly improve work efficiency.

The circulation system is made of anti-corrosion materials, and has the functions of anti-rust, anti-corrosion and low-temperature liquid pollution and sterilization.

The flow of the circulating pump can be adjusted or customized, which greatly meets the actual needs of different customers;

The microcomputer has a timer shutdown function. After enabling this function, the timer off time can be set arbitrarily from 0 to 100 hours.

The product is easy to use and can work continuously.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265