

Product Name : Medical Liquid Nitrogen Cylinder	Product Code : ASTRALAB-SUPPLYA12001
---	--



Description :
Medical Liquid Nitrogen Cylinder
Technical Specification :
The liquid nitrogen cylinder is designed for small liquid nitrogen storage and short-distance transportation, with a small neck and low evaporation loss rate.
The liquid nitrogen container is made of high-strength, special corrosion-resistant aluminum alloy, which is light in weight and easy to carry; the container is well-made, excellent in insulation, integrally formed, small welding,

firm and reliable, low loss, and the vacuum quality is more than 5 years



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265