

**Product Name :**  
Medical Liquid Nitrogen Dewar

**Product Code :**  
ASTRALAB-SUPPLYA12002



**Description :**

Medical Liquid Nitrogen Dewar

**Technical Specification :**

10 Litre portable liquid nitrogen tank is mainly used for making liquid nitrogen ice cream. It's easy to carry, simple operation, easy to transport and other advantages.

Liquid nitrogen tanks can generally be divided into storage tanks and transportation tanks. The storage tank is mainly used for the static storage of indoor liquid nitrogen, and it is not suitable for long-distance transportation

under working conditions; the transportation tank has a special shock-proof design to meet the transportation conditions.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265