

Product Name :
Medical Liquid Nitrogen Dewar

Product Code :
ASTRALAB-SUPPLYA12002



Description :

Medical Liquid Nitrogen Dewar

Technical Specification :

10 Litre portable liquid nitrogen tank is mainly used for making liquid nitrogen ice cream. It's easy to carry, simple operation, easy to transport and other advantages.

Liquid nitrogen tanks can generally be divided into storage tanks and transportation tanks. The storage tank is mainly used for the static storage of indoor liquid nitrogen, and it is not suitable for long-distance transportation

under working conditions; the transportation tank has a special shock-proof design to meet the transportation conditions.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265