

Product Name : Medical Liquid Nitrogen Container and Tank	Product Code : ASTRALAB-SUPPLYA12003
---	--



Description :
Medical Liquid Nitrogen Container and Tank
Technical Specification :
The liquid nitrogen container is mainly used for the static storage of indoor liquid nitrogen, and is not suitable for long-distance transportation under working conditions; the liquid nitrogen transportation tank has a special shock-proof design in order to meet the transportation conditions



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265