

Product Name :
Medical High Pressure Vessel

Product Code :
ASTRALAB-SUPPLYA12005



Description :

Medical High Pressure Vessel

Technical Specification :

The high-pressure container has good insulation, aluminum alloy is durable, and liquid nitrogen can hold for a long time. The container can store long biological specimens, preserve animals, semen, and biomedical vaccines.

With 6 filter tanks for biological sample storage, it is only suitable for storing liquid nitrogen and not suitable for transportation.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265