

Product Name : Medical Ice Lined Refrigerator	Product Code : ASTRALAB-SUPPLYA13002
---	--



Description :
Medical Ice Lined Refrigerator
Technical Specification :
4-digit LED high-brightness digital display, allowing users to set the temperature within the range of 2~8?, and the temperature display accuracy can reach 0.1?
Equipped with environmentally friendly CFC refrigerant R290 and compressors provided by internationally renowned brands, the refrigerator has the characteristics of fast cooling and low noise.

Built-in door handle can prevent damage during transportation and save space

There are 4 universal wheels and 2 leveling feet, which can be easily moved and fixed



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265