

**Product Name :**  
Medical Combined Refrigerator and Freezer

**Product Code :**  
ASTRALAB-SUPPLYA13004



**Description :**

Medical Combined Refrigerator and Freezer

**Technical Specification :**

Combination refrigerator with highest 2°C~8°C lowest -10°C~-26/-40°C

Electronic temperature control can achieve high precision and constant temperature

Separate cooling method, the upper fan is directly cooled, and the lower fan is directly cooled

Equipped with 2 compressors to control the upper refrigerating compartment and the lower freezing compartment respectively

Digital temperature display can accurately control the temperature and clearly monitor the operating status

Equipped with stainless steel internal material and three-layer stainless steel partition

Equipped with sheet metal refrigerated drawers and acrylic panels

Independent door locks and independent external padlocks between rooms to ensure safe storage of samples

Electric heating and LOW-E, double glass, good anti-condensation effect

The LED lighting in the combined refrigerator-freezer cabinet is easy to operate

Standard built-in USB data logger for temperature data recording

Three-layer ABS drawers in the refrigerator compartment;

There are 2 adjusting feet at the bottom and 4 universal wheels at the front and rear, which are easy to move and place



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265