

Product Name : Medical Ethylene Oxide Sterilizer	Product Code : ASTRALAB-SUPPLYA14001
--	--



Description :
Medical Ethylene Oxide Sterilizer
Technical Specification :
Ethylene oxide sterilization equipment is a special equipment for low-temperature fumigation and sterilization of items enclosed in the sterilization chamber using pure ethylene oxide at a certain temperature, pressure, humidity and certain action time.
The principle of ethylene oxide sterilization is to release a very specific alkylation effect from ethylene oxide and

microorganisms, thereby inhibiting the normal reaction set of microorganisms and causing them to suffer from metabolic disorders and death. Ethylene oxide gas is a very active chemical disinfectant, which not only achieves broad-spectrum sterilization, but also complete sterilization.

Plastic products, books, chemical fiber products, wood products, ceramic products, hysteroscope, laparoscope, laryngoscope, probe, rigid endoscope, tube endoscope, ophthalmic lens, fiber optic cable, intracranial pressure sensor, cold therapy probe , Prostatectomy device, esophagus dilator, electrocautery equipment, defibrillation electrode, laser handpiece, metal



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265