

<b>Product Name :</b> Medical Horizontal Autoclave	<b>Product Code :</b> ASTRALAB-SUPPLYA15002
-------------------------------------------------------	------------------------------------------------



<b>Description :</b>
Medical Horizontal Autoclave
<b>Technical Specification :</b>
The horizontal autoclave uses gravity exchange to more completely release the cold air from the chamber to ensure reliable sterilization. During the disinfection process, the control system will automatically adjust the steam inlet and outlet according to the temperature in the box. This equipment is ideal for clinics, scientific research institutions and other organizations to sterilize surgical instruments, fiberglass fabrics and culture media.

Automatically control the sterilization process, easy to operate.

With drying function, suitable for drying medical dressings.

With overheating, overpressure automatic protection function.

Do not open the door until the pressure in the chamber drops to 0.027MPa. If door cannot be closed well, the machine cannot be started.

When the internal pressure exceeds 0.24MPa, the safety valve will automatically open and the steam in the water tank will be discharged.

Automatically cut off the power when there is water shortage, and send out an alarm at the same time.

The chamber of the sterilizer is made of stainless steel.

With printer.



## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265