

Product Name : Medical Front Loading Autoclave	Product Code : ASTRALAB-SUPPLYA15004
--	--



Description :
Medical Front Loading Autoclave
Technical Specification :
The Front loading autoclave uses gravity exchange to more completely release the cold air from the chamber to ensure reliable sterilization. During the disinfection process, the control system will automatically adjust the steam inlet and outlet according to the temperature in the box. This equipment is ideal for clinics, scientific research institutions and other organizations to sterilize surgical instruments, fiberglass fabrics and culture media.

Automatically controlled the sterilizing process, easy to operate.

With drying function, suitable for the medical dressing drying.

With overtemperature, overpressure autoprotect.

The door can not be open until the pressure in the chamber is reduced to 0.027MPa. And the unit can not be started on if the door doesn't close well.

The safety valve will be open automatically when the inner pressure over 0.24MPa and the steam would be exhausted the water tank.

Automatically cut off the power if the waterlack happened, and alarm at the meanwhile.

Fully stainless steel structure.

With printer.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265