

Product Name :
Medical Table Top Autoclave

Product Code :
ASTRALAB-SUPPLYA15005



Description :

Medical Table Top Autoclave

Technical Specification :

The places where table top autoclave are used are commonly found in: general hospitals, stomatology, operating rooms, ophthalmology, dental clinics and beauty hospitals.

Quick disinfection within 4~6 minutes.

The screen is a touch digital display.

Has a water cut protection system.

Exhaust h 3 fixed disinfection program cycles, which can be adjusted by the user.

The whole process of adding water, heating, disinfecting, and drying steam discharge is automatically controlled.

Steam internal circulation system: no steam emission, and the environment to be disinfected will be clean and dry.

Automatic safe cooling of air.

Door safety lock system.

With four stainless steel sterilization plates.

The chamber of the sterilizer is made of stainless steel.

After disinfection, it will automatically shut down and beep.

With drying function.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265