

Product Name :
Medical Industrial Ultrasonic Cleaner

Product Code :
ASTRALAB-SUPPLYA16005



Description :

Medical Industrial Ultrasonic Cleaner

Technical Specification :

The industrial ultrasonic cleaner is specially designed for chemical reaction, chemical synthesis, biopharmaceutical synthesis, biological material degradation, recombination and Chinese medicine extraction. Develop special instruments for temperature control and cleaning of precision instrument parts. It has ultrasonic, imported closed fluorine-free compressor, 304 stainless steel heater, 304 stainless steel tank, microcomputer control, large screen LCD, drain valve, basket, stainless steel. Cover and other components.

The instrument is easy to use, ultrasonic power, ultrasonic time and reaction temperature can be set arbitrarily. The working parameters, power-off memory, and temperature accuracy of $\pm 2^{\circ}\text{C}$ are beyond the reach of ordinary ultrasound, and can meet the experimental requirements of the experimenter.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265