

Product Name :
Astra Automatic Spectrophotometer

Product Code :
ASTRALAB-SUPPLYA17002



Description :

Astra Automatic Spectrophotometer

Technical Specification :

The spectrophotometer has become a routine instrument in modern molecular biology laboratories. It is often used for nucleic acid, protein quantification and quantification of bacterial growth concentration.

Automatic detection, accurate measurement

Built in high-speed CPU. 2G memory, fast response, no need to wait

Integration of tablet computer, user-friendly interface, measurement as long as a touch!

All-in-one, working independently, one device meets all requirements

Ultimate industrial design

Precision work, more stable and reliable, simple operation.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265