

Product Name : Astra Nanodrop Spectrophotometer	Product Code : ASTRALAB-SUPPLYA17003
---	--



Description :
Astra Nanodrop Spectrophotometer
Technical Specification :
I Micro sample detection: the sample on the base only needs 0.5-2 ? L, saving precious samples.
I Automatic detection: the automatic detection function can be realized by opening and closing the upper arm.
I Wide detection range: the detection concentration range of the base is 2-27500ng / ? L

I Easy to operate: built in 7 inch high resolution color capacitive touch screen, friendly operation interface

I The best choice for nucleic acid detection: including all wavelengths required for nucleic acid concentration and purity detection.

I Multiple detection modes: it has two detection modes of micro and cuvette, which can easily deal with various types of samples.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265