

**Product Name :**  
Astra Nanodrop Spectrophotometer

**Product Code :**  
ASTRALAB-SUPPLYA17003



**Description :**

Astra Nanodrop Spectrophotometer

**Technical Specification :**

I Micro sample detection: the sample on the base only needs 0.5-2  $\mu$  L, saving precious samples.

I Automatic detection: the automatic detection function can be realized by opening and closing the upper arm.

I Wide detection range: the detection concentration range of the base is 2-27500ng /  $\mu$  L

I Easy to operate: built in 7 inch high resolution color capacitive touch screen, friendly operation interface

I The best choice for nucleic acid detection: including all wavelengths required for nucleic acid concentration and purity detection.

I Multiple detection modes: it has two detection modes of micro and cuvette, which can easily deal with various types of samples.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265