

<b>Product Name :</b> Astra Handheld Particle Counters	<b>Product Code :</b> ASTRALAB-SUPPLYA17008
---	--



<b>Description :</b>
Astra Handheld Particle Counters
<b>Technical Specification :</b>
The instrument adopts button-type one-key measurement, and high-quality spectral data can be obtained in a few seconds.
The instrument has a cable-free handheld system design, which can provide high-precision, full-band NIR spectrum measurement. The instrument adopts button-type one-key measurement, and high-quality spectral

data can be obtained in a few seconds. Optimized the spectral data collection process, more portable, simpler to operate, waterproof (but should be moisture-proof as much as possible), built-in whiteboard reference and wavelength calibration, to ensure that the instrument can be used when it is turned on. It is an ideal tool for establishing a spectral library to support image classification management and material quantitative analysis.

Full band (350-2500nm)

Ultra-portable handheld design, light and strong

Fast and accurate spectrum acquisition

No cables, no need to connect to a computer to use

Integrated computer and LCD display

Built-in whiteboard reference and wavelength calibration, ready to use

With GPS and recording function, real-time recording of measurement location and voice notes

Equipped with sample holder for scanning thin and light samples



## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265