

Product Name :
Astra Benchtop XRF Analyzer

Product Code :
ASTRALAB-SUPPLYA17011



Description :

Astra Benchtop XRF Analyzer

Technical Specification :

The benchtop xrf analyzer inherits more than 20 years of XRF research and development experience, and introduces digital multi-channel technology to make the detection limit lower, higher stability, wider application range, and performance comparable to desktop computers; small and portable size makes the detection easier ,easier.

More convenient operation

Light weight, small size, ergonomic handle design, equipped with a special instrument cover, easier to grasp, more convenient to use in the field. To

360 ° rotatable 5-inch high-definition screen, supports multi-point control, can be displayed clearly in any light. To

Sealed integrated design, with waterproof and dustproof functions, can be used continuously in harsh environments. To

There is no need to prepare samples and can directly measure the surface of the object to be tested. The instrument can be hand-held for quick testing, and it can also be used to test samples for a longer period of time.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265