

**Product Name :**  
Astra Visible Spectrophotometer

**Product Code :**  
ASTRALAB-SUPPLYA17012



**Description :**

Astra Visible Spectrophotometer

**Technical Specification :**

The sample room of the instrument can place cuvettes of various specifications between 5-50mm, which is convenient to use.

It is suitable for the basic laboratory of environmental protection major in universities and colleges as a teaching instrument.

The unique optical system, high-performance holographic grating and receiver are designed to ensure that the instrument has excellent performance indicators;

Single machine operation, linear regression equation measurement, etc.;

Apply the latest microcomputer processing technology to make the operation more convenient;

Manual wavelength adjustment;

Working method, T/A/C/F.

Standard configuration:

1 host

Power cord 1

1 cm 4-slot manual cuvette holder 1

1 cm standard glass cuvette 1 set (4 pcs)

Host instruction manual 1

Dust cover 1



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265