

Product Name : Astra Pipet Fillers	Product Code : ASTRALAB-SUPPLYA18004
--	--



Description :
Astra Pipet Fillers
Technical Specification :
Compatible with most of the plastic and glass pipettes from 0.1 -100mL
Six speeds of selection for aspiration and dispensing different liquids
Large LCD display showing low battery warning and speed settings
Enables single-handed operation with minimum effort
Light and ergonomic design providing easy usability
High capacity Li-ion battery enable long runtime of operation
Powerful pump fills a 25mL pipette in <5 seconds

0.45µm replaceable hydrophobic filter

Rechargeable during use



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265