

Product Name : Diagun	Product Code : ASTRALAB-SUPPLYA20001
--------------------------	---



Description :
Diagun
Technical Specification :
6F20 Card
12V Voltage

Display LCD with touch screen

Ambient Temperature

Relative Humidity

Performance Speed:

The whole unit (permanently connected) means the communication time between Smart box and main unit is reduced.

Appearance:

The has a new designed and cover made from high quality materials to stand in a new market.

Connector: The has a new connector to replace all CAN Bus and OBD Connector called OBDII 16E connector.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265