

Product Name :
Auto AC Recycling Machine

Product Code :
ASTRALAB-SUPPLYA20027



Description :

Auto AC Recycling Machine

Technical Specification :

Easy to read LCD display

Fast recovery

Complete recovery up to 95%.

Complete vacuum up to 28 hg

Gas charging at required quantity

Fully automatic operation

Sequential automatic operation

Fully programmable for customize operation

Auto cut on recovery, vacuum, recharging, oil separation, air purge & oil charge

Separate menu pages for simple operation

All operations can be operate individually in manual mode

Heavy duty in built pressure compressor for instant leak detection

All microprocessor controlled for error free operation

Smallest machine for space saving

Easy to maintain

Dimension: H35xW19xD24 Inch

Oil Separator: Stainless Steel

Charging Hoses: 6 Feet

Weight: 85 kg

Chaises: Powder Coated, Sheet Metal

Power: 220 V AC 50 hurts

Compressor: 1/5 sealed hermetic

Display: LCD Illuminated

Controls: Membrane Key Pad

Storage Tank Capacity: 10 kg

Pressure Gauges: 3 No. Up to 500 psi



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265