

Product Name :
Auto Steam Vacuum Cleaner

Product Code :
ASTRALAB-SUPPLYA20031



Description :

Auto Steam Vacuum Cleaner

Technical Specification :

Air Aromatization: Yes

Voltage (V): 230

Weight (kg): 20kg

Automation Grade: Automatic

Boiler Capacity (L): 2

Pressure (Bar): 5

Boiler Pressure (bar): 5

Tank Capacity (l): 2

Boiler Temperature (deg. Celsius): 190

Frequency (Hz): 50

Temperature (deg. Celsius): 2150

Cleaning Type: Steam

Boiler Power (W): 3kw

Number Of Wheels: 4

Cable Length (m): 5

Power Source: Electric

Color: Silver

Temperature: 148 Degree C

Power: 1400W

steam Vacuum cleaner is a professional steam generator and generates a pressure of 5 bar of steam in continuous fashion. this steam cleaner generates steam at 150deg centigrade. the machine also has a feature of Vacuum with heavy duty Vacuum motor which can be used with water curtain as a filter where all the dust can be collected. so you can actually steam and Vacuum at the same time with the help of different adapters provided with machine.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265