

Product Name :
PIC16F Training Lab

Product Code :
ASTRALAB-SUPPLYA21009



Description :

PIC16F Training Lab

Technical Specification :

This controller contains most of the powerful functions in modern MCUs nowadays. Moreover, it can be used for automation, motor control, device measurement, and mechanical controls...etc. PIC16F Training Lab which uses Microchip's PIC16F887 is a 8-bit microcontroller to implement various I/O control experiments. PIC16F is popular and well-known due to its economic cost, wide applicability, high accessibility and reliable stability. Together with our user friendly experiment manual, end users can learn the control of PIC MCU more conveniently and

efficiently .

PIC16F is popular and well-known due to its economic cost, wide applicability, high accessibility and reliable stability. The integration of these devices enable end users to create and implement different kind of control experiments. Contains several peripheral devices from basic LED to advanced capacitive sensing module.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265