

Product Name :
Differential Active Probe

Product Code :
ASTRALAB-SUPPLYA21042



Description :

Differential Active Probe

Technical Specification :

Bandwidth: 6 GHz.

Attenuation ratio: 2:1 or 10:1 (switchable).

Input Impedance:

Differential input R: 200 k Ω .

Single-ended input R: 100 k Ω .

Common mode input R: 50 k Ω .

Input C: 700 fF (with browser).

Range:

Dynamic range: ± 5 V or 10 Vpp (at 10:1), ± 1 V or 2 Vpp (at 2:1).

DC offset range: ± 15 V.

Maximum input voltage: 30 V DC + peak AC.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265