

**Product Name :**  
Training Bench

**Product Code :**  
ASTRALAB-SUPPLYA21045



**Description :**

Training Bench

**Technical Specification :**

A useful to create water testing dynamic. can move openly, proceed around your desk, and never having to get up from the chair. Dynamic off-ice benches can be accessed by deciding upon the ones that possess wheels. In addition, take to also to decide on those which can be corrected. This element may help it become easier that you get the most acceptable elevation of the test benches system instrumentation that is suitable for work desk. Selecting a backrest that can adapt the form of your spine is also important. Long term bench springs do not accommodate the shape of the spinal column. So, try to pick a bench with a backrest which may adhere to the pure shape of your spinal column. Like a outcome, once you lean, your spine will curl and also make you tired, sore, and also potentially expertise spinal deformities.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265