

<b>Product Name :</b> Tensile Testing Equipments	<b>Product Code :</b> ASTRALAB-SUPPLYA21064
---	--



<b>Description :</b>
Tensile Testing Equipments
<b>Technical Specification :</b>
Tensile tests up to 20 kN on specimens made of different metals, to find material characteristics such as upper and lower yield strengths, tensile strength and overall extension. Tests of Modulus for the specimen material.
Operating Conditions:

Operating environment: Laboratory.

Operating relative humidity range: 80% at temperatures < 31°C decreasing linearly to 50% at 40°C.

Storage temperature range: °C to +55°C (when packed for transport).

Operating temperature range: +5°C to +40°C.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265