

Product Name :
Planetary Mono Mill Machine

Product Code :
ASTRALAB-SUPPLYA21066



Description :

Planetary Mono Mill Machine

Technical Specification :

The cover is fitted with an inlet and outlet valve with quick-action vent. Additional lock-system – If the grinding bowl will be filled in a glove box, the additional lock-system must be used for the transport of the full grinding bowl.

Using a special lid for the grinding bowl, material can also be ground in an inert atmosphere. The Planetary

Mono Mill Pulverisette 6 is designed for fine grinding of laboratory samples or solids, dry or in suspension, down to colloidal fineness.

The Planetary Ball Mill Pulverisette 6 is compatible with the Gas and Pressure Temperature Measuring System – GTM. There is a radio transmitter in the lid and a receiver transfers the data to a computer. Continual monitoring of gas pressure and temperature enable thermal effects, physical and chemical reactions (pressure increase or decrease) to be monitored “in situ” in the grinding bowl.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265