


Product Name :
Creep Testing Machine

Product Code :
ASTRALAB-SUPPLYA21084



Description :

Creep Testing Machine

Technical Specification :

Flat specimens made of lead and PE. Tensile stress range 5...25N/mm². Table unit for the investigation of the creep behaviour of different materials. Specimen cross-section 2x5mm, measured length 25mm. Specimen cross-section 2x5mm, measured length 25mm. Position measurement gauge 0...10mm, resolution 0.01mm.

Experiments above or below room temperature possible with conditioning box

Tensile stress range: 5...25N/mm².

Position Measurement.

Measured length: 25mm.

Resolution 0.01mm.

Range: 0...10mm.

Material: lead, PE.

Cross-section: 2x5mm.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265