

**Product Name :**  
Compression Testing Machine Electrically Operated

**Product Code :**  
ASTRALAB-SUPPLYA21086



**Description :**

Compression Testing Machine Electrically Operated

**Technical Specification :**

The Digital Compression Testing Machine has been designed to meet the need for a simple, economic and reliable means to test concrete for its compressive strength. Relays are incorporated so that the motor does not restart on its own after a power breakdown. The Digital Indicator incorporates a 4-Digits display calibrated in Kilo Newton (KN) and

preset to maximum load capacity fitted with micro switches to switch-off the motor when the load approaches the maximum capacity of the gauge, avoiding any over loading. The electrically operated pumping units are provided with a control knob to adjust the pace rate which can be effectively controlled by an experienced operator during the course of testing, by observing the progress of the Digital Indicator readings.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265