

Product Name :
Compression Testing Machine Electrically Operated

Product Code :
ASTRALAB-SUPPLYA21086



Description :

Compression Testing Machine Electrically Operated

Technical Specification :

The Digital Compression Testing Machine has been designed to meet the need for a simple, economic and reliable means to test concrete for its compressive strength. Relays are incorporated so that the motor does not restart on its own after a power breakdown. The Digital Indicator incorporates a 4-Digits display calibrated in Kilo Newton (KN) and

preset to maximum load capacity fitted with micro switches to switch-off the motor when the load approaches the maximum capacity of the gauge, avoiding any over loading. The electrically operated pumping units are provided with a control knob to adjust the pace rate which can be effectively controlled by an experienced operator during the course of testing, by observing the progress of the Digital Indicator readings.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265